Analyzing the Relationship Between the Acidity and the Conductivity of a Solution

Tarik Onalan

23 February 2015

1 Introduction

1.1 Purpose

Understand how pH of a solution affects conductivity measured in μS .

1.2 Hypothesis

Conductivity will decrease as pH approaches 7, as the ions that would otherwise facilitate the conduction of electricity (H^+/OH^-) would be neutralized.

1.3 Variables

Independent Variable

• pH

Dependent Variable

 \bullet Conductivity $[\mu S]$

Controlled Variables

- Molarity of HCl
- Molarity of NaOH
- Volume of HCl (50 mL)
- Drop rate of NaOH ($\sim 1 \, \mathrm{s}^{-1}$)
- Spinning speed of stir bar

2 Materials

- 150 mL 0.10 m HCl
- $\bullet~150\,\mathrm{mL}~0.10\,\mathrm{M}$ NaOH
- $1 \cdot 100 \,\mathrm{mL}$ buret+stand

- $1 \cdot 250 \,\mathrm{mL}$ beaker
- $\bullet~1~\cdot~50\,\mathrm{mL}$ graduated cylinder
- 1 · stir plate+bar
- $1 \cdot \log ging device$
- $\bullet~1\cdot \mathrm{pH}$ probe
- 1 · conductivity probe
- $\bullet~1~\cdot~\mathrm{USB}$ flash disk
- Phenolphthalein

3 Procedure

- 1. Turn on logger
- 2. Plug in pH/conductivity probes and USB disk
- 3. Set up buret with stand
- 4. Put stir plate below buret
- 5. Fill graduated cylinder with 50 mL HCl
- 6. Fill buret with 50 mL NaOH
- 7. Transfer HCl to beaker
- 8. Put 2-3 drops of phenolphthalein in the beaker
- 9. Calibrate probes
- 10. Put beaker on top of stir plate
- 11. Put the stir bar in the beaker
- 12. Set the stir bar spinning at $\frac{1}{4}$ speed
- 13. Open the burst so there is $\sim 1 \text{ drop s}^{-1}$
- 14. Start logger
- 15. Let experiment run for 50 seconds
- 16. Stop logger
- 17. Export data to USB disk
- 18. Release any remaining NaOH into HCl solution
- 19. Turn off stir plate
- 20. Remove stir bar
- 21. Clean beaker
- 22. Repeat 5-21 as necessary for data collection

4 Data




Table 1: Trial 1			
Time	pH Conductivity		
\mathbf{S}		μS	
0.00	1.01±0.20	19 352.00	±800.
2.00	1.02 ± 0.20	19374.00	$\pm 800.$
4.00	1.04 ± 0.20	19323.00	$\pm 800.$
6.00	1.07 ± 0.20	19293.00	$\pm 800.$
8.00	1.11 ± 0.20	19272.00	$\pm 800.$
10.00	1.08 ± 0.20	19301.00	$\pm 800.$
12.00	1.14 ± 0.20	19264.00	$\pm 800.$
14.00	1.17 ± 0.20	19293.00	$\pm 800.$
16.00	1.17 ± 0.20	19257.00	$\pm 800.$
18.00	1.19 ± 0.20	19257.00	$\pm 800.$
20.00	1.21 ± 0.20	19257.00	$\pm 800.$
22.00	1.24 ± 0.20	19228.00	$\pm 800.$
24.00	1.34 ± 0.20	19220.00	$\pm 800.$
26.00	1.27 ± 0.20	19176.00	$\pm 800.$
28.00	1.35 ± 0.20	19191.00	$\pm 800.$
30.00	1.35 ± 0.20	19139.00	$\pm 800.$
32.00	$1.42 {\pm} 0.20$	19125.00	$\pm 800.$
34.00	1.39 ± 0.20	19088.00	$\pm 800.$
36.00	1.37 ± 0.20	19030.00	$\pm 800.$
38.00	$1.42 {\pm} 0.20$	18985.00	$\pm 800.$
40.00	$1.47 {\pm} 0.20$	18949.00	$\pm 800.$
42.00	$1.45 {\pm} 0.20$	18934.00	$\pm 800.$
44.00	1.50 ± 0.20	18905.00	$\pm 800.$
46.00	1.49 ± 0.20	18846.00	$\pm 800.$
48.00	1.49 ± 0.20	18839.00	$\pm 800.$
50.00	$1.57 {\pm} 0.20$	18780.00	$\pm 800.$

	Table 2: Trial 2		
Time	рН	Conductivi	ity
\mathbf{S}		μS	
0.00	$0.82 {\pm} 0.20$	19 646.00	±800.
2.00	0.79 ± 0.20	19669.00	$\pm 800.$
4.00	$0.83 {\pm} 0.20$	19646.00	$\pm 800.$
6.00	$0.84 {\pm} 0.20$	19661.00	$\pm 800.$
8.00	$0.86 {\pm} 0.20$	19639.00	$\pm 800.$
10.00	$0.87 {\pm} 0.20$	19632.00	$\pm 800.$
12.00	0.90 ± 0.20	19609.00	$\pm 800.$
14.00	$0.94 {\pm} 0.20$	19579.00	$\pm 800.$
16.00	$0.98 {\pm} 0.20$	19579.00	$\pm 800.$
18.00	$0.95 {\pm} 0.20$	19579.00	$\pm 800.$
20.00	0.97 ± 0.20	19541.00	$\pm 800.$
22.00	1.01 ± 0.20	19527.00	$\pm 800.$
24.00	0.98 ± 0.20	19504.00	$\pm 800.$
26.00	1.04 ± 0.20	19474.00	$\pm 800.$
28.00	1.03 ± 0.20	19451.00	$\pm 800.$
30.00	1.08 ± 0.20	19421.00	$\pm 800.$
32.00	1.05 ± 0.20	19414.00	$\pm 800.$
34.00	1.06 ± 0.20	19369.00	$\pm 800.$
36.00	1.13 ± 0.20	19331.00	$\pm 800.$
38.00	1.11 ± 0.20	19286.00	$\pm 800.$
40.00	1.11 ± 0.20	19264.00	$\pm 800.$
42.00	1.13 ± 0.20	19226.00	$\pm 800.$
44.00	1.17 ± 0.20	19226.00	$\pm 800.$
46.00	1.19 ± 0.20	19174.00	$\pm 800.$
48.00	$1.21 {\pm} 0.20$	19144.00	$\pm 800.$
50.00	$1.22 {\pm} 0.20$	19099.00	$\pm 800.$

Table 3: Trial 3			
Time	рН	Conductivi	ity
\mathbf{S}		μS	
0.00	0.40 ± 0.20	19 046.00	±800.
2.00	$0.50 {\pm} 0.20$	19053.00	$\pm 800.$
4.00	0.50 ± 0.20	19024.00	$\pm 800.$
6.00	$0.56 {\pm} 0.20$	18988.00	$\pm 800.$
8.00	$0.58 {\pm} 0.20$	18988.00	$\pm 800.$
10.00	$0.61 {\pm} 0.20$	18959.00	$\pm 800.$
12.00	$0.63 {\pm} 0.20$	18930.00	$\pm 800.$
14.00	$0.71 {\pm} 0.20$	18872.00	$\pm 800.$
16.00	$0.76 {\pm} 0.20$	18814.00	$\pm 800.$
18.00	$0.81 {\pm} 0.20$	18778.00	$\pm 800.$
20.00	$0.85 {\pm} 0.20$	18727.00	$\pm 800.$
22.00	0.89 ± 0.20	18662.00	$\pm 800.$
24.00	$0.96 {\pm} 0.20$	18604.00	$\pm 800.$
26.00	$1.05 {\pm} 0.20$	18517.00	$\pm 800.$
28.00	$1.05 {\pm} 0.20$	18481.00	$\pm 800.$
30.00	1.09 ± 0.20	18401.00	$\pm 800.$
32.00	1.14 ± 0.20	18285.00	$\pm 800.$
34.00	1.17 ± 0.20	18213.00	$\pm 800.$
36.00	$1.25 {\pm} 0.20$	18155.00	$\pm 800.$
38.00	$1.26 {\pm} 0.20$	18076.00	$\pm 800.$
40.00	1.27 ± 0.20	17981.00	$\pm 800.$
42.00	1.31 ± 0.20	17865.00	$\pm 800.$
44.00	1.29 ± 0.20	17779.00	$\pm 800.$
46.00	$1.31 {\pm} 0.20$	17685.00	$\pm 800.$
48.00	1.20 ± 0.20	17612.00	$\pm 800.$
50.00	$1.22 {\pm} 0.20$	17482.00	$\pm 800.$

Table 4: Average

\overline{Time}	pH Conductivity		
\mathbf{s}		μS	
0.00	0.74 ± 0.20	19348.00	$\pm 800.$
2.00	0.77 ± 0.20	19365.33	$\pm 800.$
4.00	0.79 ± 0.20	19331.00	$\pm 800.$
6.00	$0.82 {\pm} 0.20$	19314.00	$\pm 800.$
8.00	$0.85 {\pm} 0.20$	19299.67	$\pm 800.$
10.00	$0.85 {\pm} 0.20$	19297.33	$\pm 800.$
12.00	0.89 ± 0.20	19267.67	$\pm 800.$
14.00	$0.94 {\pm} 0.20$	19248.00	$\pm 800.$
16.00	0.97 ± 0.20	19216.67	$\pm 800.$
18.00	$0.98 {\pm} 0.20$	19204.67	$\pm 800.$
20.00	1.01 ± 0.20	19175.00	$\pm 800.$
22.00	1.05 ± 0.20	19139.00	$\pm 800.$
24.00	1.09 ± 0.20	19109.33	$\pm 800.$
26.00	1.12 ± 0.20	19055.67	$\pm 800.$
28.00	1.14 ± 0.20	19041.00	$\pm 800.$
30.00	1.17 ± 0.20	18987.00	$\pm 800.$
32.00	1.20 ± 0.20	18941.33	$\pm 800.$
34.00	1.21 ± 0.20	18890.00	$\pm 800.$
36.00	$1.25 {\pm} 0.20$	18838.67	$\pm 800.$
38.00	$1.26 {\pm} 0.20$	18782.33	$\pm 800.$
40.00	1.28 ± 0.20	18731.33	$\pm 800.$
42.00	1.30 ± 0.20	18675.00	$\pm 800.$
44.00	$1.32 {\pm} 0.20$	18636.67	$\pm 800.$
46.00	1.33 ± 0.20	18568.33	$\pm 800.$
48.00	1.30 ± 0.20	18531.67	$\pm 800.$
50.00	$1.34 {\pm} 0.20$	18453.67	$\pm 800.$

Calculations 5

5.1 Average

$$\langle \mu S(pH_i) \rangle$$
 (1)

$$\frac{\sum_{i} \mu S(pH_i)}{i} \tag{2}$$

$$\frac{\mu S(pH_1) + \mu S(pH_2) + \dots + \mu S(pH_{i-1}) + \mu S(pH_i)}{i}$$
(3)

5.2 Pearson Correlation Coefficient

The Pearson Correlation Coefficient $r_{x,y}$ is defined as

$$r_{x,y} = \frac{\sum_{i=1}^{n} cov(x,y)}{\sigma_x \sigma_y} \tag{4}$$

where cov(x,y) is the covariance and σ_x/σ_y are the standard deviations of sets x and y. This can be expanded to

$$r_{x,y} = \frac{\sum_{i=1}^{n} (x_i - \overline{x})(y_i - \overline{y})}{\sum_{i=1}^{n} (x_i - \overline{x})^2 \sum_{i=1}^{n} (y_i - \overline{y})^2}$$
(5)

where $\overline{x}/\overline{y}$ are the averages of sets x and y. $r_{x,y}$ was then calculated by substituting sets x and y with the pH and conductivity sets, respectively.

5.3 Error and Minimization

Error was calculated using sum squared error, which is defined as follows:
$$E = \sum_{i=1}^{n} (y_i - f(x_i))^2 \tag{6}$$

While a program carried out the minimization of the function, the basic premise of minimization is to "follow the gradient," so to speak. The program follows the gradient

$$E' = \frac{\partial}{\partial x} \left(\sum_{i=1}^{n} (y_i - f(x_i))^2 \right) \tag{7}$$

similar to Euler's method for solving differential equations. Minimizing the function is simply following the "negative" slope to the local minima of the function.

5.4 Line of Best Fit

While creating my best fit line, I took into account the axis units. pH is a logarithmic scale, and µS is a linear scale. From there, I hypothesized that my line of best fit would be logarithmic in nature.

$$f(x) = kln(px) + h (8)$$

However, this did not produce lower error than a second order polynomial,

$$q(x) = ax^2 + bx + c (9)$$

which was surprising. The root cause seemed to come from the fact that the logarithmic curve, by nature, "flattens out" much more rapidly than a polynomial curve. To compensate, I concatenated the two functions:

$$(f \circ g)(x) = kln(p(ax^2 + bx + c)) + h \tag{10}$$

$$(f \circ g)(x) = k\ln(ax^2 + bx + c) + h \tag{11}$$

This approach led to a much lower error than either function alone, as shown below:

Error			
	$(f \circ g)(x)$	f(x)	g(x)
	27423.187126	286923.387218	46056.2186406

6 Conclusion

My hypothesis, that the conductivity would decrease as pH approached 7, was partially correct. I was only able to test with an acidic starting point, meaning that I only observed the pH approaching 7 from $0 \le \text{pH} < 7$. However, I can conclude that as the pH approaches 7 from $0 \le \text{pH} < 7$, the conductivity decreases. Simply from a mathematical standpoint, the Pearson correlation coefficient—the measure of the "correlatedness" of any two datasets—of the pHand conductivity returns -0.95, which is an almost exact (negative) linear relationship. This is evident by looking at the graph, which shows that as the pH increased, the conductivity decreased. When the pH was ~ 0.74 (starting pH), the conductivity was $19\,348.00\,\text{pS}$. When the pH increased to ~ 1.34 , however, the conductivity decreased to $18\,453.70\,\text{µS}$. Of course, there are many more values in between, but the general trend is similar to the one described.

The biggest difficulty I had in my lab was keeping the probes calibrated. In my lab, I was using 0.10 M HCl and NaOH, which have pH values of 1 and 13, respectively. However, the pH probe in particular would not hold its calibration, and in my second and third trials, the starting pH of the HCl is incorrectly reported as 0.82 and 0.40. Even so, this does not invalidate my conclusion, because there is still an obvious downward trend in the conductivity-pH data.

In the future, I could compensate for probe calibration issues by running more trials, so that I could average away systematic error. In this lab, I only had enough time to run 3 trials, but, in the future, I could run more trials to compensate for sensor drift.